


National Distributor

Komet Medical.
















EZblade® – Premium saw blades for TKA.



This premium knee saw blade made by Komet Medical is provided with a structured surface, reduced edges and teeth with specially ground clearance.



Explore our Catalogue here.

Widths	Cutting Thicknesses	Connections	Compatible Connections/Hub Styles
 <p>13 mm</p>	 <p>1,00 mm</p>	 <p>Aesculap 3Ti</p>	<p>Acculan 3Ti (GA673, GA668), Acculan 4 (GA331, GA336)</p>
 <p>19 mm</p>	 <p>1,19 mm</p>	 <p>DePuySynthes/ Zimmer Biomet</p>	<p>DePuy Synthes Compact Air Drive II + Power Drive (511.800, 511.801), Air Power Line (511.610), Battery Power Line (530.610), Battery Power Line II (530.710), TRS (05.001.224, 05.001.223, 05.001.240)</p> <p>Zimmer Biomet Zimmer Universal Power System (89-8509-450-60)</p>
 <p>22 mm*</p>	 <p>1,27 mm</p>  <p>1,37 mm</p>	 <p>ConMed/ Linvatec/Hall</p>	<p>M Power (PRO6300), M Power 2 (PRO6350M), PowerPro Max (PRO5300M), PowerPro + PowerPro Electric II (PRO6125), PowerPro Pneumatic (PRO6175), Hall Serie 4 (5067-002), Hall Serie 3 (5044-002), Versipower Plus (5071-002), Hall 50 (PRO7300B)</p>
 <p>25 mm</p>	 <p>1,47 mm</p>	 <p>MicroAire</p>	<p>Series 7000</p>
		 <p>Stryker</p>	<p>System EHD (296-98), System 2000 (2108), System 4 (4108), System 5 (4208), System 6 (6208), System 7 (7208), System 8 (8208)</p>
		 <p>DeSoutter</p>	<p>MDX Multi Drive (KX 600), MDX Multi Drive (KDX-600), MPX Multi Drive (SX 500), MPX Multi Drive (DPX-676), MCX Multi Drive (DCX 500), Orthodrive (DBK-700) Not compatible with DBK-701</p>

*available in parallel and flared

Contact:

Allegra Orthopaedics Limited

Level 8, 18-20 Orion Road Lane Cove NSW 2066

T +612 9119 9200 Toll Free 1800 644 370

E sales@allegraorthopaedics.com

W www.allegraorthopaedics.com

